Table: Face-to-face Communication – some options for people with a range of needs:

Communication Devices	Mainstream Apple /Android Devices with Speech Generating Apps	Specialist Devices with optional keyguard and/or personalised settings	Specialist Devices with Eye Gaze Cameras	Specialist Devices with Head Mouse Cameras	Specialist Device used with Switch/Joystick /Mouse
Who they are for	Good for people who can use their hands or a stylus to make selections / use a keyboard	Good For people who can use their hands or a stylus, who need specialist symbol vocabulary (usually long term AAC users)	Good for people who are unable to use their hands or a stylus to make selections/ use a keyboard	Good for people who are unable to use their hands or a stylus to make selections / use a keyboard	Good for people who are unable to use their hands or a stylus to make selections / use a keyboard
How they work	Touch to type or Touch or to select words/ phrases/symbols/ Pictures	Touch to select words/phrases/ symbols/pictures or Touch to type	 Eye gaze devices enable hands-free communication. People use their eye movements, tracked by a camera, to navigate round a device and make selections. Options are available for people who use text or symbols, 	 A head mouse enables handsfree communication. Slight movements of the head mouse (attached as a small dot to the forehead or glasses) are tracked on screen by a camera. This enables navigation round a device to make selections. 	 A switch/joystick of mouse can be used to scan and select text or symbols on a communication device. Switches, joysticks and mice come in all shapes and sizes and can be operated by any part of the body from your head to your feet. Options are available for people who use text or

			or a combination of both.	Options are available for people who use text or symbols, or a combination of both.	symbols, or a combination of both
Advantages	 Other functions a available of the iPad/ Android, s as email a social med A small number of apps can used with banked voices and messages 	 Other functions can be added such as email/social media. In Kent & Medway these are added by the Environmental Control Service. 	 Can be used hands-free. Other functions can be added such as email/ social media. In Kent & Medway these are added by the Environmental Control Service. Some models are compatible with banked voices and/or messages. 	 Can be used hands-free, when head and neck movement is comfortable and reasonably strong. Other functions can be added such as email/social media. In Kent & Medway these are added by the Environmental Control Service Some models are compatible with banked voices and/or messages. 	 Can be used if hands/arms are limited or uncontrolled. Other functions can be added such as email/ social media. In Kent & Medway these are added by the Environmental Control Service. Some models are compatible with banked voices and/or messages.